



RAMPART ARMERY

Rampart's Armery radial system enables easy and consistent left and right radial access while attenuating radiation and protecting the team.

Radiation Attenuating

Closes gap to minimize radiation scatter by up to 65%1

Adjustable

The left radial accessory provides easy access and aids patient comfort during left radial procedures

Simple and Standardized Left and Right Access

The Armery system has been customdesigned to fit with the Rampart Bunker Table Shield, providing interventional staff with a >1mm barrier of protection at the table level.



Innovative Supportive Work Bench

Runway supports wires, catheters, and other procedural equipment

Lead-free







